

Fuzion® Cleaner Disinfectant



Bleach efficacy that's easy on surfaces

Clorox Healthcare® Fuzion® Cleaner Disinfectant combines the proven power of bleach with excellent aesthetics for broad use throughout the facility.

- ▶ **Kills C. diff in 2 Minutes:** Hypochlorous acid solution is a strong oxidizer and highly effective surface disinfectant.
- ▶ **Broad Surface Compatibility:** Safe for use on stainless steel, clear plastics, finished wood and many other common surfaces.
- ▶ **Low Residue:** Leaves surfaces looking clean, not cloudy.
- ▶ **Low Odor:** Mild odor dissipates in minutes.
- ▶ **Easy to Use:** No PPE required.

| Fuzion® at a Glance | |
|--|----------|
| <i>C. difficile</i> contact time | 2 min |
| TB contact time | 1 min |
| Bacteria, virus and fungi contact time | 1 min |
| Kills Norovirus? | Yes |
| Signal word | Caution† |
| HMIS rating | 0-0-0-X |
| PPE required | None |

*Use as directed on hard, nonporous surfaces. Follow special instructions for cleaning prior to disinfection.

†Toxicity Signal word on Product Label. Fuzion is EPA-registered and exempt from GHS requirements.

Clorox Healthcare® Fuzion® Cleaner Disinfectant Spray


| Effective Against: |
|--|
| Spores - 2 min kill time |
| <i>Clostridioides difficile</i> spores |
| Bacteria - 1 min kill time |
| <i>Acinetobacter baumannii</i> |
| <i>Bordetella pertussis</i> |
| <i>Bordetella bronchiseptica</i> |
| <i>Burkholderia cepacia</i> |
| <i>Campylobacter jejuni</i> |
| <i>Enterobacter cloacae</i> (NDM-1) |
| <i>Enterococcus faecalis</i> |
| <i>Enterococcus faecalis</i> (VRE) (Vancomycin-resistant) |
| <i>Enterococcus faecium</i> (MDR) (Multidrug-resistant) |
| <i>Enterococcus hirae</i> |
| <i>Escherichia coli</i> O157:H7 (<i>E. coli</i>) |
| <i>Klebsiella aerogenes</i> |
| <i>Klebsiella pneumoniae</i> (CRE) (Carbapenem-resistant) |
| <i>Klebsiella pneumoniae</i> (NDM-1) |
| <i>Legionella pneumophila</i> |
| <i>Listeria monocytogenes</i> |
| <i>Proteus mirabilis</i> |
| <i>Pseudomonas aeruginosa</i> |
| <i>Salmonella enterica</i> |
| <i>Serratia marcescens</i> |
| <i>Shigella dysenteriae</i> |
| <i>Staphylococcus aureus</i> |
| <i>Staphylococcus aureus</i> (MRSA) (Methicillin-resistant) |
| <i>Staphylococcus aureus</i> (MDR) (Multidrug-resistant) |
| <i>Staphylococcus aureus</i> (VISA) (Vancomycin-intermediate) |
| <i>Stenotrophomonas maltophilia</i> |
| <i>Streptococcus pneumoniae</i> |
| <i>Streptococcus pneumoniae</i> (MDR) (Multidrug-resistant) |
| <i>Streptococcus pyogenes</i> |

| Viruses - 1 min kill time |
|--|
| Adenovirus Type 2 |
| Adenovirus Type 5 |
| Avian Influenza A Virus (H5N1) |
| Canine Adenovirus |
| Canine Parvovirus |
| Coxsackie Type B3 Virus |
| Enterovirus D68 |
| Feline Panleukopenia Virus |
| Feline Rhinotracheitis Virus |
| Hepatitis A Virus |
| Hepatitis B Virus [†] |
| Hepatitis C Virus* |
| Herpes Simplex Virus Type 1 (HSV-1) |
| Human Coronavirus |
| Human Immunodeficiency Virus Type 1 (HIV-1) |
| Influenza A Virus |
| Influenza B Virus |
| Measles Virus |
| Middle East Respiratory Syndrome Coronavirus (MERS-CoV) |
| Murine Parvovirus |
| Norovirus [‡] |
| Poliovirus Type 1 (Chat) |
| Respiratory Syncytial Virus (RSV) |
| Rhinovirus |
| Rotavirus |
| SARS-Associated Coronavirus |
| SARS-CoV-2 (cause of COVID-19) |
| Fungi - 1 min kill time |
| <i>Candida albicans</i> |
| <i>Candida glabrata</i> |
| <i>Microsporium canis</i> |
| <i>Trichophyton interdigitale</i> (Athlete's Foot Fungus) |
| TB - 1 min kill time |
| <i>Mycobacterium bovis</i> (TB) |
| TOTAL: 66 CLAIMS |
| Use as directed on hard, nonporous surfaces. |

Excellent surface compatibility

Using Fuzion® eliminates the need to make tradeoffs between efficacy and surface compatibility. Fuzion® is a pH-neutral formula that won't cause damage to common surfaces. It can be safely used on stainless steel, chrome, countertops, sinks, glazed tiles — even clear plastics and glass.

In external laboratory testing, Fuzion® left significantly less residue on a clear acrylic surface compared to a quat/alcohol and a peracetic acid/hydrogen peroxide disinfectant:



Fuzion® enables facilities to use bleach in more places, providing strong protection against the spread of pathogens on surfaces.

* Obtained via Duck Hepatitis B Virus surrogate.

† Obtained via Bovine Viral Diarrhea Virus surrogate.

‡ Obtained via Feline Calicivirus surrogate.

NI-57842

Product Information

Clorox Healthcare® Fuzion® Cleaner Disinfectant

EPA Reg. No. 67619-30
Spray Bottle, 32 oz., 9/case
Item No. 31478
Case UPC No.
10044600314782



HEALTHCARE®

FROM CLOROX PRO

For product resources and implementation tools,
contact your Clorox sales representative or
Call: 1-800-234-7700
Visit: www.CloroxHealthcare.com

© 2022 Clorox Professional Products Company,
1221 Broadway, Oakland, CA 94612.